



AP BLOT

AUTOMATED BLOT PROCESSOR



USER MANUAL

This manual contains DAS reserved information. All rights reserved. Unauthorized copying of this manual or parts of it is prohibited.



INDEX

SECTION 1 INTRODUCTION	7
1.1. DESCRIPTION	7
1.2. WARRANTY	7
1.2.1. Requesting service.....	7
1.2.2. Ordering spare parts.....	7
1.3. CE CONFORMITY CERTIFICATE	8
SECTION 2 - GENERAL SAFETY.....	9
2.1. GENERAL SAFETY WARNING SYMBOLS	9
2.2. INSTRUMENT WORKING CONDITIONS.....	9
2.3. INTENDED USE OF THE INSTRUMENT.....	9
SECTION 3 - GENERAL INFORMATION	11
3.1. INSTRUMENT DESCRIPTION: PURPOSE AND FEATURES	11
3.2. MAIN FEATURES	11
3.3. TECHNICAL SPECIFICATIONS.....	11
SECTION 4 - INSTALLATION AND START-UP	13
4.1. UNPACKING.....	13
4.2. POSITIONING THE INSTRUMENT	15
4.3. POWER CONNECTION	15
4.4. INSTALLING WASH AND WASTE CONTAINERS.....	15
4.5. WORK AREA.....	16
4.6. PC CONNECTION	16
4.7. INSTALLING APEX SOFTWARE ON PC.....	17
4.7.1. Minimum PC requirements	17
4.7.2. Installing APEX software.....	17
4.7.3. Installing XCALIB software	17
4.7.4. Installing USB driver	18
4.7.4.1. Windows XP	18
4.7.4.2. Windows 2000.....	19
4.7.5. Installing camera uEye 3.20 driver.....	22
SECTION 5 – INSTRUCTIONS FOR USE.....	25
5.1. GENERAL INDICATIONS.....	25
5.1.1. Introduction.....	25
5.1.2. Turning instrument on.....	26
5.1.3. Running APEX programme.....	26
5.1.4. LICENCES.....	26
5.1.4.1. Installing licences.....	26
5.1.5. Main Menu	30
5.2. SESSION MENU	32
5.2.1. NEW SESSION.....	32
5.2.2. Work Area Display.....	49
5.2.2.1. Running a session.....	55
5.2.2.2. Displaying Session Results.....	57
5.2.3. RUN LAST SESSION.....	59
5.2.4. LOAD FROM HOST	59
5.2.5. EXIT.....	59
5.3. MENU DATABASE	60
5.3.1. ASSAY.....	60
5.3.2. LAYOUTS.....	60
5.3.2.1. New layout/Set-up	61
5.3.2.2. Changing Layouts	62
5.3.2.3. Deleting a set-up/layout	62
5.3.2.4. Cloning a layout.....	63
5.3.2.5. Renaming layouts	63
5.3.2.1 Configuration protection.....	64
5.3.2.2 Report Print Out.....	65
5.3.2.6. Importing layouts	65
5.3.2.7. Exporting layouts	65
5.3.3. RESULT ARCHIVES.....	66
5.4. MENU SETTINGS	68
5.4.1. LANGUAGE.....	68
5.4.2. INSTRUMENT.....	68
5.4.3. USERS.....	69
5.4.3.1. NEW USER.....	70
5.4.3.2. MODIFY.....	71



5.4.3.3.	DELETE	71
5.4.3.4.	CHANGE PASSWORD.....	71
5.4.4.	CHANGE USER.....	72
5.4.5.	CHANGE PASSWORD	73
5.5.	LABORATORY DATA.....	74
5.5.1.	DOCTOR DATA.....	75
5.5.2.	IMPORT SESSION DATA.....	76
5.5.3.	EXPORT SESSIONS DATA.....	76
5.5.4.	SETTING COMMUNICATION PORT.....	77
5.5.5.	COVER BLOCK.....	78
5.5.6.	OPTIONS.....	79
5.6.	TOOLS MENU	87
5.6.1.	MOVEMENT TESTS.....	87
5.6.2.	PRIME TEST.....	88
5.6.3.	LIGHT ON.....	89
5.7.	MENU ?	90
5.7.1.	Guide.....	90
5.7.2.	LICENSES.....	90
5.7.3.	INSTRUMENT INFORMATION.....	92
5.7.4.	SOFTWARE COMPONENTS.....	93
5.7.5.	ABOUT.....	94
SECTION 6 – SPECIALIST ACCESS		95
6.1.	SPECIALIST and ADVANCED SPECIALIST.....	95
6.2.	SETTING UP A NEW ASSAY.....	95
6.2.1.	Phases Settings.....	98
6.2.2.	Changing an Assay.....	105
6.3.	Copying an existing assay	106
6.4.	Printing out an assay	106
6.5.	Importing a protocol	106
6.6.	Exporting an assay	108
SECTION 7 – MAINTENANCE		109
7.1.	MAINTENANCE CHECKS & PRECAUTIONS.....	109
7.1.1.	Cleaning.....	109
7.1.2.	Instrument testing.....	109
7.1.3.	Inspections and checks.....	109
7.2.	REPAIRS	110
SECTION 8 - PACKING, TRANSPORTATION, STORAGE AND PUTTING INSTRUMENT “OUT OF SERVICE” .		111
8.1.	GENERAL INFORMATION.....	111
8.2.	PUTTING INSTRUMENT “OUT OF SERVICE”/DECOMMISSIONING	111
8.2.1.	Temporary Storage.....	111
8.3.	PACKING THE INSTRUMENT	111
8.4.	TRANSPORTATION.....	112
8.5.	STORING THE INSTRUMENT	112
8.6.	DISPOSAL PROCEDURES.....	112
SECTION 9 – APPENDIX		113
APPENDIX A.....		113
A.1.	INSTRUMENT DECONTAMINATION	113
A.1.1.	DECONTAMINATION PROCEDURE.....	113
APPENDIX B.....		115

SECTION 1 INTRODUCTION

1.1. DESCRIPTION

This manual provides Users with all the necessary instructions and maintenance recommendations for a safe and suitable use of the instrument.

Manual content:

- **Section 1** - warranty information, Technical Service and the CE conformity declaration
- **section 2** - general safety-warnings
- **section 3** - general information, instrument performance and technical data
- **section 4** - installation and start up
- **section 5** - operating instructions for IFA process
- **section 6** - instructions for IFA specialist
- **section 7** – instructions for regular maintenance by User(s)
- **section 8** - how to put instrument out of service, packaging, transport, storage and decommissioning
- **section 9** - appendices

This manual is considered an integral part of the instrument; it should always be available for User/Maintenance Technician consultation.

For accurate installation, usage and maintenance please read all the instructions carefully before use.

In particular carefully read Section 1 “GENERAL SAFETY WARNINGS”, describing suitable operating procedures to avoid risk of harm or damage to instrument and personnel.

In case of problems or instrument breakdown, refer to Authorised Technical Service and Repairs only.

1.2. WARRANTY

Each provided instrument is completely tested, and guaranteed for twelve months from delivery.

The warranty is for all mechanical and electrical parts. It is granted only for a proper installation, use and maintenance in compliance to the given instructions contained in this manual.

Due to a misuse of the instrument the above mentioned warranty does not include the travel and labour-hour expenses as well as all the accommodation expenses.

The warranty does not include any responsibility for direct or indirect to people and material damages, caused by improper use and maintenance of the instrument. All parts that are subject to deterioration, because of their specific use, are excluded from the warranty.

1.2.1. Requesting service

When REQUESTING TECHNICAL SERVICE please refer to this manual and indicate any data displayed on instrument identification label (model, registration number and year of manufacture).

Please Note: only qualified Technicians with suitable equipment and original spare parts are entitled to repair the instrument; ordinary maintenance should be carried out by Users as described in this manual.

1.2.2. Ordering spare parts

Parts that are subject to deterioration or that are defective and need to be replaced have to be requested with the following information:

- Customer purchase order number
- Name/type and/or version of the instrument
- instrument code number
- part code number
- part description
- requested quantity
- name and company address for delivery of ordered original spare parts

Replacing parts with ORIGINAL SPARE PARTS guarantees effective and longer lasting instrument life.

1.3. CE CONFORMITY CERTIFICATE



Declaration of conformity in compliance with 98/79/EC In-Vitro Diagnostic medical devices Directive (IVD)

DAS srl, manufacturer of the Automated Blot Processor “**AP-Blot**”, In-Vitro Diagnostic Medical Devices risk class, declares that this instrument complies with the European Directive 98/79/EC (IVD).

The present conformity declaration, as stated in Annex III of the 98/79/EC Directive, is based upon the following requirements:

- Quality Management System (ISO 9001:2000): Working Procedures and Instructions for
 - Design process control;
 - Manufacturing process control;
 - Tests, Checks and inspections.
- System documentation:
 - Technical documents;
 - Risk analysis;
 - User manual;
 - Service manual.

The instrument fully conforms to the following EU directives, which include latest modifications and comply with national regulations and legislation.

- **Directive 2006/95/EC (ex 73/23/CE) regarding low voltage safety** and related standards:
 - EN 61010-1:2001 - Safety Precautions for electrical, measuring, control and laboratory equipment
- **Directive 2004/108/CE (ex 89/336/EC) - directive regarding electromagnetic compatibility** and related standards:
 - CEI EN 61326-1 (for electrical, measuring, control and laboratory equipment)
 - EN 55011 (Emission RF)
 - EN 61000-3-2 (Harmonics)
 - EN 61000-3-3 (Flickers)

Palombara Sabina, October 27th, 2008

Mario Orsini
Managing Director

SECTION 2 - GENERAL SAFETY

2.1. GENERAL SAFETY WARNING SYMBOLS

In this manual the following symbols stand for danger or warnings:



General DANGER symbol. Indicates that not following relative instructions and warnings can be a serious safety risk.



Indicates ELECTRICAL VOLTAGE, which could cause death upon contact. Covers with this symbol can only be removed and replaced by qualified people and only after electrical power disconnection.



Indicates that instrument is used with reagents and other DANGEROUS CHEMICAL SUBSTANCES, corrosive, irritating or noxious, which could damage health.



Indicates that the instrument handles potentially infectious samples (e.g. body fluids such as urine) which could cause INFECTION/CONTAMINATION. Always observe the general safety standard rules for the presence of these biological substances.



Indicates that not following correct instructions could cause damage to the instrument and/or its proper functioning.



Indicates important information concerning the instrument or a section of the document that is important to read with particular attention.

2.2. INSTRUMENT WORKING CONDITIONS

The instrument is intended to be used in the following working conditions:

- as In Vitro Diagnostic (IVD) medical device as specified in the technical data
- with chemical reagents and accessories supplied and/or declared to be compatible with the instrument
- at a specific temperature and humidity level according to the details given in this manual
- not to be used and powered in a potentially explosive or fire hazardous environment



The instrument should only be used as described in this manual.
Any other use is regarded as improper

2.3. INTENDED USE OF THE INSTRUMENT

The instrument has to be used for the expected purpose and in perfect technical conditions, by qualified personnel, in working conditions and following maintenance operations described in this manual.

This manual contains instructions for qualified operators.

- **Qualified Users** have to ensure that environmental conditions are suitable and that installation, use and maintenance are correct and in accordance with general safety rules, as well as following those specific precautions described in this manual.
- No User is qualified or entitled to repair the instrument.
- **Only Qualified and specifically trained Technicians** are entitled to service and repair the instrument with original spare parts.



Modifications of the instrument are not allowed. The user is liable for any improper modifications and for the consequences deriving from them.

For extraordinary maintenance, ask for licensed service centers. The maintenance will be carried out by **Specialized Technicians** that will be able to repair the instrument using original spare parts.



SECTION 3 - GENERAL INFORMATION

3.1. INSTRUMENT DESCRIPTION: PURPOSE AND FEATURES

The AP-Blot instrument dispenses, rinses and reads blots on immuno-enzymatic analysis strips. It can also dispense, incubate and wash organic samples (pre-diluted when required) and specific reagents. Up to 48 strips can be processed.

3.2. MAIN FEATURES

Aim	Automatic analysis of Dot-Blot and Western Blot tests, capable of processing up to 10 different protocols simultaneously
Samples	Maximum capacity: 48 samples in 13x100mm tubes. Programmable Automatic Predilution for each sample directly carried out on trays
Reagents	Modular Tray for controls, conjugates, chromogens, stop solutions, common reagents (e.g.: sample dilution solutions)
Reaction Areas	6 racks for 8 different positions: 48 strip total capacity
Dispensation	Three axis arm (X, Y, Z ₁ , Z ₂) with two independent probes, one with liquid level sensor
Washing	<i>Probe washing:</i> external immersion wash in-well; internal wash with Buffer solution. <i>Strip:</i> shaking plane washing system
Image acquisition	High resolution colour camera
External liquid bottles	A probe washing WASH bottle container, a WASH BUFFER container for blot washing and a WASTE bottle container. Each bottle has a level sensor. Each container is inside the instrument
Immuno-enzymatic management software	PC software is Windows XP and 2000 compatible capable of: – Setting up work lists – Automatically executing all session phases (dispensation, incubation, washing and readings) – Creating final reports and result archive

3.3. TECHNICAL SPECIFICATIONS

Sample rack	48 sample tubes (Ø 13mm)
Reagent rack	6 x 100ml phials and 8 control positions
Liquid Containers	2 bottles for wash buffer and waste. Both with liquid level sensors
Strip Holders	6 slides (8 strips per slide) to process up to 48 strips on shaking plane
Dispensation Probe Unit	1 dispensation probe and 1 aspiration probe
Strip Washing	Liquid dispensation and aspiration via two instrument probes
Image acquisition	High Resolution Colour Camera (3.3 Megapixels) mounted above Probes
Dispensation	Up to 2500 µl with a 1 µl resolution
Software	Windows 2000 and XP compatible to set up work list, reports and result archive. Loading work lists and transmit results via Host Computer interface (LIS)
Minimum PC Requirements	– 200 MB free hard disk space – CD-Rom reader – RAM 512 MB – Pentium 4 – Standard video card:VGA 1024x768 16 million colours – Free USB 2.0 port
Power Requirements	230/115 VAC, 50-60Hz, 150W
Optional	Bar Code Reader (BCR) for samples with position ID
Measurements (in cms)	79 x 61 x 65 h
Weight	78 Kg

SECTION 4 - INSTALLATION AND START-UP




When installing and setting up the instrument, carefully follow safety warnings and general rules described in **Section 1**.


4.1. UNPACKING

The instrument is packed with a wooden box. To unpack the instrument follow the instructions described below:


- a) Remove the 2 plastic covers (1-2) from the instrument right side and the other ones from the left side. Screw in the 4 handles on both sides (3-4) for manual instrument placement. Once in position remove the handles and put back in place the plastic covers.
- b) Remove the 4 screws that hold the instrument on to the wood base (5-6)



- c) Remove screw and bracket (7), screw in the same place using the screw removed before



- d) Remove screws (8,9) and related brackets
- e) Remove the Z block holder (10) as indicated
- f) Remove the protection cap from the camera (11)
- g) Remove the blocking tools (12) from the oscillating plane



If the instrument was already in use check for the **DECONTAMINATION DECLARATION** (see Appendix A of this manual).



KEEP THE EMPTY BOX! This box is designed and manufactured specifically for this instrument only and can be used for future transportation.

4.2. POSITIONING THE INSTRUMENT

The instrument has to be placed on a level bench (table or equivalent work surface) and in an environment where the temperature ranges between 15 ° C and 35 ° C, humidity levels remain under 80% and the instrument is protected from direct sun light.

4.3. POWER CONNECTION

Once the instrument has been placed in position proceed with plugging it into the correct power source 230VAC or 115VAC (see label) using the cables supplied with instrument.



Warning: check that electrical power source corresponds to the source indicated on the label stuck on the side of the instrument!



Warning: check that fuses correspond to those indicated on the side of the instrument!

230 VAC n° 2 x 5A fast fuses

115 VAC n° 2 x 10A fast fuses

4.4. INSTALLING WASH AND WASTE CONTAINERS

Figure 1 shows how wash/waste containers are connected up to the instrument.




fig. 1 - Hydraulic circuit

- R** Level sensors connector
- T1** **Wash:** probe wash solution bottle
- T2** **Buffer:** strip wash solution bottle
- T3** **Waste:** waste liquid bottle

PROCEDURE

Before switching on instrument:

- fill T1 and T2 bottles with solution according to method requirements
- connect empty T3 (waste) bottle container
- connect tubing to appropriate containers as indicated in figure 1
- check tubings are correctly connected and in good condition
- plug level sensor lead to instrument socket (connector R)



Note: at the end of a session or before performing maintenance procedures: run wash session for instrument.

4.5. WORK AREA




fig. 2- AP-Blot Work Area

4.6. PC CONNECTION

The instrument is equipped with a USB port for PC connection. To install the USB driver please see par. 4.7.4.




fig. 3- Electrical socket/ 2.0 USB port



4.7. INSTALLING APEX SOFTWARE ON PC

4.7.1. Minimum PC requirements

- Microsoft Windows 2000 or XP
- 200 MB free hard disk space
- CD-Rom reader
- 512 MB RAM
- Pentium 4
- Standard video card VGA 1024x768 16 million colours
- 2.0 USB port

4.7.2. Installing APEX software

Insert "APEX software" CD into computer. APEX software is self-installing. If this automatic installation process does not begin, it can be launched manually by clicking on **setup.exe** file which appears on the display.

Choose programme display language (English or *Italiano*). NB: *this choice is only valid for the installation.*

Click "Next" until installation procedure window appears:

- Click on "Complete" to choose default file - C:\Programmes\APEXSuite\APEX;
- Click on "Personalise" to change default file destination path.

After installing APEX software continue to install USB driver as described later (see para 4.7.4).

At the end of the installation launch APEX software by double clicking on APEX icon which appears on the desktop.

NB: User can access all APEX suite programmes from the "APEX" group from main menu (by selecting Start>all programmes>APEX Suite).

4.7.3. Installing XCALIB software

NB: XCALIB software is the AP-Blot Service software

Insert "XCALIB software" CD into computer. XCALIB software is self-installing.

If this process does not automatically start then click on the **setup.exe** file from the CD displayed.

Choose language after installation has begun (English or *Italiano*). NB: *this choice is only for the installation procedure.*

Click on "Next" until installation choice window is displayed:

- Click on "Complete" to install default file C:\Programmes\APEXSuite\XCALIB;
- Click on "Personalise" to change installation file path destination.

When software is installed continue by installing USB driver as described later (see para 4.7.4).

At the end of installation procedure, launch XCALIB software by clicking on XCALIB icon which appears on the desktop as part of the "APEX Suite" group. Or From the main menu programme bar (select "start>all programmes>APEX Suite)

4.7.4. Installing USB driver

The instrument communicates with PC through a convertor USB port (RS232-USB).

2 driver installation procedures are required: the first to communicate with the electronic part of the instrument and the second to emulate (=imitate) a standard PC serial port.

When these drivers have been correctly installed a new serial port has been added to the PC (e.g.: COM5).

4.7.4.1. Windows XP

Connect instrument to the PC USB. The “**New Hardware Guided Installation**” procedure will automatically start up with the following display:




fig. 4

Select “No, not yet” in answer to question: “Connect to Windows Update to search for new software?” then click on “Next”

Select “Install from list or specific path” then click on “Next”.




fig. 5

Click “**BROWSE**” to find file containing USB driver, then click **OK** and continue.

If Windows XP is configured to warn of the installation of a non-WHQL certified driver a “Warning! Message” will appear: click on “Continue anyway”.

The Windows programme displays a note at the end of the installation - click on “**END**” to complete.

The same procedure now automatically begins to install a serial port emulator – follow instructions as above.

When both drivers have been installed, open Control Panel select System then click on Hard Ware file then **Hardware Resources**.

The new USB/Serial port should appear as an additional serial port called “USB Serial Port” (see below for example – called “COM7” in fig. 6).




fig. 6

Now APEX or XCALIB software can be launched and the appropriate communication port can be selected.

4.7.4.2. Windows 2000

Connect instrument to PC USB. The “New Hardware Guided Installation” procedure will automatically start up. Click on “Next” to continue.



Fig. 7

Select "Search for appropriate Driver"
Continue by clicking on "Next".



fig. 8

Select only "Specify path" and click on "Next".



fig. 9

Click on “BROWSE” and specify USB driver path file in display which appears below:




Fig. 10

When file with USB driver has been opened, select **FTDIBUS.INF** and then click “Next” to continue installation.

A Windows display announces successful installation and the USB is ready.

A similar process begins for the serial port emulator driver. Follow similar procedure.

When both drivers have been installed open Control Panel, select System and then Hardware file and open the “**Hardware Resources**”.

The “new” USB/Serial Port should appear as an additional communication port: “USB Serial Port” (e.g.: COM5 as shown below Fig.12).




Fig. 11

4.7.5. Installing camera uEye 3.20 driver

4.7.5.1. Windows XP

To install uEye 3.20 camera driver proceed as follows:

- 1) Disconnect instrument from PC 2.0 USB
- 2) From driver CD click on display which appears: "CAMERA DRIVER uEye_3.20" then select file "uEyeBatchInstall_320.exe".




fig. 12

- 3) Select "Install the uEye USB driver" option by typing in '2'. Press Enter




fig. 13

- 4) When installation is successfully completed a confirmation message will appear: "uEye driver installed correctly". To continue press Enter.




fig. 14

- 5) Now connect APBLOT instrument to the PC 2.0 USB port.
- 6) When new hardware is recognised "Connect to Windows Update to search for new software?" appears and now select "No, not yet" - then click on "Next" to continue.



fig. 15

- 7) Select "Automatic software installation" and click "Next".




fig. 16

- 8) If a display shows incompatibility with Windows XP click on "Continue" and installation continues:




fig. 17

- 9) The following display will appear when installation has been successful:



fig. 18

- 10) Click on End
11) NB: sometimes camera driver installation software requests procedure to be repeated from point (6)
12) Now launch APEX or XCALIB software and choose appropriate communication port.

SECTION 5 – INSTRUCTIONS FOR USE

5.1. GENERAL INDICATIONS



The instrument is used with chemical reagents and other DANGEROUS CHEMICAL SUBSTANCES (corrosive, irritant and harmful) which can damage health. When this label is present particular care and attention must be taken to follow all manufacturer recommendations.



The instrument handles potentially infected samples (urine and human serum etc.) therefore INFECTIONS or CONTAMINATION might occur.

Take all proper precautions in the presence of such biological substances.

Use protective clothing, gloves and goggles.

In particular:

- Wear protective clothing, goggles and gloves
- At the end of a work session, check condition and carefully clean work panel and accessories with a cloth soaked in a 0.5% sodium hypochlorite solution
- Any waste material produced during the sessions and cleaning process can be dangerous – with potential risk of contamination for personnel and other matter.
- Dispose of such waste material in accordance with local laws and regulations

5.1.1. Introduction

The following table shows the main way this instrument functions.




fig. 19– General functions

5.1.2. Turning instrument on

Press the on/off button on the right hand side of instrument. There are three lights on the front of the machine:

POWER	<u>Green light</u>	Indicates that the instrument is switched ON
RUN	<u>Yellow light</u>	Indicates a programme is running
ALARM	<u>Red light</u>	This shows that at least one liquid error has been found: the tip has not found any liquid in one of the containers (sample tube or reagent container) on the work area (Note 1)

Please note: The position where the liquid error has been found will be displayed on the PC screen.

5.1.3. Running APEX programme

The Apex programme icon appears on the desktop. Double click on it to launch programme.



Apex

Fig. 20 – Programme icon

5.1.4. LICENCES

5.1.4.1. Installing licences

Double click on programme icon to begin English language programme – the following window display appears asking if any new licenses are to be loaded (Fig. 29 below):




Fig. 21 – License request

Click **No** to exit this part of the programme

Click **Yes** to display Licenses window which lists possible licenses to be loaded as shown below Fig.22




fig. 22- Licenses Window

Please note: The license is only valid for one instrument.

Please note: Loading a Demo License for this specific instrument **overwrites any other licenses.**

Any other licenses present (if required) should be re-loaded.

Select **Load new license** to open display and CD installation folder and then click “Licenses”.

These are contained in an ***.xml file** and can be opened by selecting Apri (Fig. 23) or by clicking **Demo mode** and selecting type of instrument to function in Demo Mode – without being connected to a PC.

Choose License for instrument and click “Apri/Open”:




Fig. 23– Selecting the license



fig. 24

Click **Deselect all** to de-select license and not load it.
Click **OK** to load license.

Select **OK** to show license in the licenses window as shown below in Fig. 33





fig. 25


When **Licenses** window is closed (exit) the following form will be displayed automatically. Insert both login and password and select Enter.



The screenshot shows a window titled 'APEX'. It contains two input fields: 'LOGIN' and 'PASSWORD'. Below these fields are two buttons: 'Cancel' and 'Enter'.

Fig. 26- Software access form

After inserting login and password the **Communication port settings** window will open to input communication port setting from set of pull-down lists Fig. 26



The screenshot shows a window titled 'Communication port settings'. It contains three pull-down menus: 'Instrument' (showing 'APPADBC.97'), 'Communication type' (showing 'SERIAL'), and 'Port'. At the bottom, there are two buttons: 'Cancel' and 'Set'.

Fig. 27 - Port settings

Select **Set** to store selected communication port. This selection is confirmed as shown below Fig. 27
Select **Cancel** to exit this part of programme and leave original setting.

Click OK. The Main Menu will appear if the selected port is correct and the instrument is on.



Fig. 28– Confirmation window

A Warning message will appear (Fig. 28) asking if communication port should be changed. If selection is incorrect, or if Port is correct, but the instrument is not turned on.




Fig. 29– Change selected port request

If the instrument is turned off, turn it on and select **No**.



The Main Menu will now appear. Click on “light” icon (fig. 30) if the port is correct the light will come on.




fig. 30- Main Menu On/Off light

If instrument is on but the light does not flash ‘on’ this indicates that software is only in DEMO mode. Change port selection.

The light will turn on as soon as the correct port is selected.

Click on Light Icon to turn light off.

5.1.5. Main Menu

The Main Menu appears as below in Fig. 39:




Fig. 31 -Main Menu

The following diagram (Fig. 32) shows the operating menu and how it works:




Fig. 32– Functions Menu



Menu Functions

Menu:	Description:
<i>Session</i>	
New Session	To set up new work session
Run last session	Runs last session with previous settings
Load from host	Import of work list from Host Computer
Exit	Leaves programme
<i>Database</i>	
Result archives	Display results of all stored sessions and results for all patients examined in each session
Assay	Manages all functions concerning assays
Layout	This layout display manages all Strip functions set-ups.
<i>Settings</i>	
Language	Sets up display language for menus, windows etc
Instrument	To choose one or more instruments
Work area	Shows the work area of the selected instrument
Users	To insert/change/delete User password for programme.
Change User	To change any software User access.
Change password	To change password of current User
Laboratory data	To input information about laboratories where sessions take place.
Doctor data	To input information about Doctor responsible for sessions.
Import session data	To import session in *.mdb file format and overwrite stored sessions.
Export session data	To export sessions in *.mdb file format and store sessions in archive.
Communication port	To set up serial port to connect AP-BLOT instrument to PC.
Cover lock	Sets safety block levels for cover
Options	Access to different functions, about each or all instruments: error trace type, session set up forms, time storage algorithm, etc...
<i>Tools</i>	
Movement tests	Tests movements of: different Axes, shaking plane, primes hydraulic circuit, Optical sensor test
Prime test	Hydraulic circuit prime test from selected prime source
Light on	Testing light
<i>?</i>	
Guide	Displays User Manual in chosen language
Licenses	Display/load/deletes instrument licences
Instrument information	Displays serial number, code, card types and firmware version of connected and working instrument
Software components	Displays information about different APEX software components.
About	Information about APEX software

5.2. SESSION MENU

5.2.1. NEW SESSION

Click on **Session** from the Main Menu and select **New Session** (Fig 33) to begin setting up new work session from a new display. This Last Session display shows the last settings set-up on work area – please note – this last session has not necessarily been run, but has been at least set-up on Lay-Out display.

There are two steps to set up a work session: the first concerns choosing one or more assays to be executed by one or more instruments.

The second step is associating one or more samples to one or more assays. A Session can be more easily identified by filling in the Session Name field.




fig. 33

The Session Settings display opens where one or more patients can be associated to one or more assays (Fig.2).




fig. 34

In the “Work list” section click  icon to add an assay to the work list (Fig.34).




fig. 35

This opens the “Assay selection” display where Assays required can be selected (Fig.4) and added to work list after clicking “Confirm”.

Repeat this procedure if more than one assay is required for the work list.




fig. 36

Please Note: If more than one instrument is connected to the same PC, the instrument specifically required for the Assay can be selected from the “**Instrument**” pulldown listing in the “**Instruments**” section of the “**Assay selection**” display.

The different sections:

- **“General info”** : shows the general information about the Assay (Name, short name etc)
- **“Repetitions”** : shows how many repeats are set-up for both samples and controls (where 1= corresponds to a single sample/control dispensation)
- **“Controls”** : shows information about Controls dispensation – which takes place either before or after the samples.

The “Repetition” and “Controls” sections have default settings specifically for Dot Blot Assays and consequently cannot be changed.

Click Confirm on Assay Selection display and the Assay is added to the Work List and displayed on the Session Settings with the short name shown in the first column of the table (Fig.37).




fig. 37

Patient (Sample) Codes can also be added to the Work List in the **“Sample code”** column.

These codes can be inserted either one by one, or in groups (using “multiple selection” option).

to **Insert Single Code** (Fig.38): type in code (point a) in the Sample Code column, then press Enter.

The cursor flashes in cell below waiting for a new code to be typed in.

A Patient Card icon appears on the left (point b) and a small box appears on the right (point c) to input sample/assay association.





fig. 38

Multiple Patient Code Selection method (Fig.39): to insert a group of Patient Codes click on **“Sample code”** heading (point a) to open Multiple insert window (point b).


Type in **“Patient name”** prefix code and also the **“Patient number”** which corresponds to how many Patient Codes should be generated at the same time.



The image shows a dialog box titled "Multiple insert" with a dark gray background. Inside the dialog, there is a section labeled "Data" which contains two input fields: "Patient name:" and "Patient number:". Below these fields are two buttons: "Cancel" and "OK". The dialog box is overlaid on a light gray grid background.

fig. 39

For instance: Type "P" as the prefix code in "Patient name" text box and then type in the number '5' in the "Patient number" box and click "OK" (Fig.40).



Multiple insert

Data


Patient name: P

Patient number: 5

Cancel OK

fig. 40

In this case - 5 Patient numbers are automatically generated starting with "P" followed by an increasing number starting from 1 (P1, P2... etc see Fig.41).



APEX - [Session setting]

Session Database Settings Tools ?

Global info

Sample

Associated: 0
Not associated: 6

Request

To be processed: 0
Suspended: 0

Assay: A1 (A1)

Work list


	Sample Code	A1 AP-Blot 99	
	paziente 1	☐	
	P1	☐	
	P2	☐	
	P3	☐	
	P4	☐	
	P5	☐	

fig. 41

This work list now displays both samples and assays.
Now define a state for each Assay request to complete this session setting.
The possibilities are:

- Not processed: represented by ☐ icon
- To be processed: represented by  icon
- Suspended: represented by  icon
- Validated: represented by  icon
- Running: represented by  icon

Click to check each box that corresponds to associating the assay with that patient sample as shown in Fig.42. A yellow check confirms the request to process this sample/assay association which is then inserted in the work list to be processed.




The screenshot shows the 'APEX - [Session setting]' window. The 'Global info' section includes 'Sample' (Associated: 3, Not associated: 3) and 'Request' (To be processed: 3, Suspended: 0). The 'Assay: A1 (A1)' section shows 'General info' (Type: Dot blot, Instrument: AP-Blot 99), 'Request' (To be processed: 3, Suspended: 0), and 'Repetition' (Sample: 1, Control: 1). The 'Work list' table has columns for Sample Code, A1 AP-Blot 99, and several empty columns for other assays. The table contains rows for 'paziente 1', 'P1', 'P2', 'P3', 'P4', and 'P5'. The 'A1 AP-Blot 99' column has checkboxes for each row, with 'P2' checked. The 'Exit' and 'Validate' buttons are at the bottom.

Sample Code	A1 AP-Blot 99														
paziente 1	☒														
P1	☒														
P2	☒														
P3	☐														
P4	☐														
P5	☐														

fig. 42

Click again on  icon to add another Assay (Fig.43) if required; this icon moves to the right when a new assay is included and short name displayed in a column.



The screenshot shows the 'APEX - [Session setting]' window with two assay columns, 'A1 AP-Blot 99' and 'A2 AP-Blot 99'. The 'Work list' table has columns for Sample Code, A1 AP-Blot 99, A2 AP-Blot 99, and an 'Add Assay' icon. The table contains rows for 'paziente 1', 'P1', 'P2', 'P3', 'P4', and 'P5'. The 'A1 AP-Blot 99' column has checkboxes for each row, with 'P2' checked. The 'A2 AP-Blot 99' column has checkboxes for each row, with 'P2' checked. The 'Add Assay' icon is at the end of each row.

Sample Code	A1 AP-Blot 99	A2 AP-Blot 99	Add Assay
paziente 1	☒	☐	
P1	☒	☐	
P2	☒	☒	
P3	☐	☐	
P4	☐	☐	
P5	☐	☐	

fig. 43

There are also two ways to associate assays and samples with the **Multiple Selection** option:

- The first concerns **associating more than one patient to the same Assay** (a multiple selection of patient codes as shown in Fig.44): click on (down) arrow at the top of the column under the Assay Short name required, to select the whole column (when selected it appears in a darker blue).




fig. 44

Click on any box in this column to select this Assay for all patient sample codes (Fig.45).





fig. 45


Please note: this selection can easily be modified by simply 'de-clicking' the box corresponding to the patient code to be excluded from this Assay  (Fig.46).



Sample Code	A1 AP-Blot 99	A2 AP-Blot 99	
paziente 1	☒	☒	
P1	☒	☒	
P2	☒	☒	
P3	☐	☒	
P4	☐	☒	
P5	☐	☒	

fig. 46

- The second Multiple Selection method concerns **associating more Assays to the same patient**: by clicking on the arrow on the right of the Patient (Sample) Code the whole Row can be selected (Fig.47).



Sample Code	A1 AP-Blot 99	A2 AP-Blot 99			
paziente 1	☒	☒			
P1	☒	☒			
P2	☒	☒			
P3	☐	☒			
P4	☐	☒			
P5	☐	☒			

fig. 47

At this point (when row is selected it appears in a darker blue) click on any box to extend the selection to the entire row (Fig.48).




fig. 48

The Top Section of the Session Settings display (under main menu toolbar) shows the total and partial request counters (Fig.49).

The **“Global info”** section shows the overall information about the Session setting.

Under **“Sample”** - the **“Associated”** and **“Not associated”** figures refer to how many samples have been associated and how many have not been associated to at least one Assay.

The **“Request”** section shows how many samples have been requested for processing, or how many have been suspended, compared to the total number of requests present for this session (which correspond to the yellow checks in the boxes).

Click on Assay Short name from Work List to display any information about the **“Assay”** section:

“General info” shows the general information about an Assay (type, etc)

“Request” shows how many Assays have been requested for processing (☒) and those which have been suspended (☐) ; this is why an Assay counter can be ‘partial’.

“Repetition” shows how many repetitions are requested for samples and controls.




fig. 49

In the Assay section (Fig.49) click on ‘Previous’ arrow and ‘Next’ arrow to display information about other Assays which have been loaded for the session.

Click  icon to display the Assay selection window (Fig.50) where the default settings of the Assays are shown.




fig. 50

The order of Patient (Sample) Codes from the work list, which are then displayed on the sample rack, can be changed. For instance if the work list appears as shown in Fig. 51 after confirming the Work Area lay-out the sample rack is displayed as shown in Fig.52.




fig. 51




fig. 52

Use the two arrows along the right side of the Session Settings display: select a Patient Sample Code and then click to move selection up or click to move selection down.

For instance Fig. 53 shows that the first row (which corresponds to “Paziente 1” code) has been selected followed by a click on the Move Down arrow (Fig.53) .




fig. 53

This moves the “Paziente 1” sample code down one row – or into second position as shown in Fig.54. Click **Validate** and the Sample Rack display from the Work Area also displays this change (Fig.55).

Sample Code	A1 AP-Blot 99	A2 AP-Blot 99	+
P1	☒	☒	
paziente 1	☒	☒	
P2	☒	☒	
P3	☐	☒	
P4	☐	☒	
P5	☒	☒	

fig. 54




fig. 55

[illegible][illegible]

fig. 57

This example shows that both Sample Assay requests for the same Patient sample code have been suspended (Fig.58).

Work list					
	Sample Code		A1 AP-Blot 99	A2 AP-Blot 99	
	paziente 1				
	P1				
	P2				
	P3				
	P4				
	P5				

fig. 58

If a work session is now run then only those requested (with a yellow check) would be processed. The next session work list would only display the one "Paziente 1" patient sample code associated with the A1 and A2 Assays and showing their suspended state (Fig.59).

APEX - [Session setting]

Session Database Settings Tools ?

Global info

Sample

Associated: 1

Not associated: 0

Request

To be processed: 0

Suspended: 2

Assay: A1 (A1)

Work list

	Sample Code	A1 AP-Blot 99	A2 AP-Blot 99	
	paziente 1			

fig. 59

[illegible]


fig. 60

Cancelling a Patient sample can only take place after this procedure has been confirmed.

Click "**Validate**" to confirm Assay selections and associations. The yellow checks become green indicating that they have been 'validated' (Fig.61).

If any request is not confirmed a Warning Message appears.

Read only or **Full session** options are now possible from this 'validated' Session Setting display.




The screenshot shows the APEX [Session setting] window. At the top, there are tabs for Session, Database, Settings, and Tools. Below the tabs, there are several status boxes: Sample (Associated: 5, Not associated: 0), Request (To be processed: 8, Suspended: 0), Assay: A1 (A1), General info (Type: Dot blot, Instrument: AP-Blot 99, N° Control: 0), Request (To be processed: 3, Suspended: 0), and Repetition (Sample: 1, Control: 1). The main area is a table titled 'Work list' with columns for Sample Code, A1 AP-Blot 99, A2 AP-Blot 99, and a series of checkboxes. The table contains five rows of data for samples P1 to P5. At the bottom, there are buttons for Exit, Read only, and Full session.

Sample Code	A1 AP-Blot 99	A2 AP-Blot 99	Checkboxes											
P1	☒	☒	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
P2	☒	☒	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
P3	☐	☒	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
P4	☐	☒	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
P5	☒	☒	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

fig. 61

Lastly, information about a patient sample code can be included in the corresponding Patient Card: click on  icon on the left of Patient (Sample) Code to display a new window (Fig.62).

The Patient-ID code is already displayed automatically (corresponding to that on the Work List) and further information about the patient can be compiled in the other fields. Click Confirm to store information.




The screenshot shows the APEX [Session setting] window with the Patient Card dialog box open for sample P1. The dialog box has fields for Patient ID (P1), Group, Name, Surname, Date of birth, Sex (M/F), and Note. There are buttons for Cancel and Confirm. The background window shows the same Work list table as in Fig. 61.

fig. 62

The following fields are obligatory: Patient-ID, name, surname (family or last name) and date of birth (Fig.63). This Patient Card can be also compiled during the running phase as the Session setting can be opened at any time and information is only stored at the end of the session. When Patient Cards have been stored they can be displayed in a list and included again for another session by clicking on  icon(Fig. 64).

Click  icon to delete card and consequently any information which might be incorrect.



PATIENT CARD - P1

Patient-ID: Group:


Name:

Surname:

Date of birth: Sex: ☒ M ☐ F

Note:

fig. 63



List of patient cards

N.	Name	Last name	Birthdate	Group	Sex
1	nome paziente	cognome paziente	01/01/2000		

fig. 64

This current session setting display can be used to set-up a new work session while the software is running another session.

New requests of sample codes associated to assays can be added now to be processed later (for instance in Fig. 65 'Nuovo Paziente' is an example of a New Patient).



5.2.2. Work Area Display

When setting up the work session is completed the **Work Area** display can be accessed (Fig. 28)
There is a total of nine racks present on the work area:

- n°1 Probe Washing basin
- n°2 Reagent Racks
- n°1 Sample Rack
- n°6 Strip plates




Fig. 66 – Work area

Probe Washing Basin

The Probe Washing Basin is on the top left of the work area layout display.
Double click on it to display condition as shown here:




Fig. 67 – Probe washing basin

Click on **Report** to display a print preview of all the information contained in the Printed Report.
Click on the images of each single rack to open a display of selected rack and a list of information about each rack phial. Different windows now show different racks.

SAMPLE RACK




Fig. 68 – Sample rack

The Sample Rack is composed of 48 different interactive positions.

Click on any position to display all the information about that particular position.

The colour of each sample position shows whether a sample is present. Red = present; Grey a sample is not present.

The information along the bottom of the window:

- ☐ SAMPLE POSITION
- ☐ TYPE
- ☐ ID
- ☐ REQUIRED VOLUME

Click on **Report** at the bottom of the window to display a print preview of all the information that can be printed out.

REAGENT RACK




Fig. 69 – Standard Reagent rack with controls




fig. 70– Additional reagent rack

There are two racks, one for reagents with positions: numbered R1 through to R9, and a control rack with positions C1 to C8 and there is a CL position for a Cleaning solution.

These positions are interactive and show different colours depending on the Session assays.

The example of the assay session is composed of the following flacons:

- ☐ PRE-TREATMENT (R2)
- ☐ SAMPLES AND CONTROLS (R1 and C1)
- ☐ CONJUGATE (R7)
- ☐ SUBSTRATE (R8)
- ☐ ADDED PHASE 1 (R3)
- ☐ ADDED PHASE 2 (R4)
- ☐ ADDED PHASE 3 (R5)
- ☐ ADDED PHASE 4 (R6)
- ☐ STOP SOLUTION (R9)

The positions for the washing solution is also shown (CL).
Eight other positions show the controls.

The colour indicates whether the control is present: green = present; grey is not present.

Click on any position in the table or representation on the display, to show all the information regarding that container
The information displayed on the left of the window:

- ☐ POSITION
- ☐ TYPE
- ☐ ID
- ☐ REQUIRED VOLUME
- ☐ ASSAY

Click on **Report** at the bottom of the window to display a print preview of all the information that can then be printed out.

STRIP PLATE




Fig. 71 – Strip Plate

There are 8 strips for each strip plate.

Information about any of the eight strips on the strip plate can be displayed simply by clicking on the selected strip.

Different colours represent either sample or control strips:

green represents a control strip, red represents a sample strip.

Information about any specific strip or position can also be displayed by clicking on the strip position (in the Data section under the Position listing).

The information displayed along the bottom of the window are:

- ☐ POSITION
- ☐ TYPE
- ☐ IDENTIFIER
- ☐ REQUIRED VOLUME
- ☐ ASSAY

Different strips can be used in all six strip plates to execute different tests at the same time, although all the strip plates are structurally the same.


Click on **Report** at the bottom of the window to display a print preview of all the information that can then be printed out.

Along the bottom of the Work Area display: the **Next** option opens the next page and the **Previous** option opens the previous page.

Kit Data, **Report** and **Info Buffer** are explained as follows:

KIT DATA

Select **Kit Data** option to display a window to input information about each assay kit fig 72.




The 'Kit data' window is a dark-themed dialog box with a title bar containing a small icon and the text 'Kit data'. It features two dropdown menus at the top: 'Assay' set to 'All the assays' and 'Kit index' set to 'All indexes'. Below these is a section titled 'Data' containing several input fields: 'ID' (text box), 'Lot' (text box), 'Producer' (text box), 'Expiry date' (calendar dropdown showing '21/05/2009'), 'Assay group' (text box with 'G1'), 'Assay' (text box with 'A1'), 'ShortName' (text box with 'A1'), and 'Assay type' (text box with 'DotBlot'). At the bottom are 'Cancel' and 'Save' buttons.

Fig. 72 – Kit Data

Kit Information can be referred to one or more assays and they can be numbered. Each kit information can include: ID number, Lot number, Producer, Expiry Date.

Cancel option closes this window without any additional information and **Save** stores all data which will also be displayed with the results (fig. 73).



The 'Results' window is a dark-themed application window with a title bar showing a small icon and the text 'Results'. It contains a 'Session' section with 'Date' (14/07/2009), 'Time' (09:26), and 'Number of assays' (3). To the right is an 'Instrument' dropdown set to 'AP-Blot (s.n. 99)'. Below the session info is an 'Assay' list on the left showing '1 A11', '2 A12', and '3 A1'. The main area on the right is titled 'GENERAL INFORMATION' and contains a table with the following data:

Session	New
Assay group	G1
Assay	A11
Execution	14/07/2009 09:26:39

Below this table is a section titled 'KIT INFORMATION' and another titled 'PROTOCOL INFORMATION'. Under 'PROTOCOL INFORMATION', there is a sub-section 'Interval interpretation' showing a horizontal scale from 10 to 60, with markers for 'Negative', 'Borderline', 'Low positive', and 'High positive'.

fig. 73– Results Information

BUFFER INFO

Click **Buffer Info** along the bottom of the Work Area display to display the Buffer Report showing which Buffer has been used to wash strips (fig. 74):




fig. 74– Buffer Info Report

REPORT

Click on **Report** along the bottom of the Work Area display to show a print preview with the following information about slides and racks:

- ☐ RACK POSITION
- ☐ TYPE
- ☐ ID
- ☐ REQUIRED VOLUME
- ☐ ASSAY




fig. 75– Report

5.2.2.1. Running a session

Click **Next** option on Work Area window (Fig. 66) to open **Running Session** window (Fig. 76).
The analysis sessions start here and results are displayed at the end of each session.




Fig. 76 – Running session

The assays used in the analysis session are listed at the top of the window. For each assay a timeline shows the progress throughout the work phase. The right side of the screen shows the assay name, the current phase of the assay, and the remaining time in hours-minutes-seconds.

The horizontal bar in the centre of the screen shows the percentage or rate of progress of each given session.

Directly underneath this bar four sections show: first any information about the **Current Session**, followed by second section to display any **Errors** that might occur during the session, the third shows information about **Thermal Regulation**, fourth shows the **Session Time** – how long before it finishes or: if completed it shows when it started and finished.

The **Unblock lid** option is to open the shield cover during a session if the safety level is high.

The **Clean list** option empties the errors list present. All errors are still printed out in the final report.

Select **Start** to run analysis session.

The **Pause** key can temporarily put the session on hold. To restart click on **Start**.

Any analysis session can be definitively abandoned by selecting the **Stop** button.

The **Back** option displays previous window, even during session.

The **Exit** option closes run session phase. It is not active during the session.

ASSAY INFO

Click on Assay name (Fig.77) to display information about set-ups as shown here:




Fig. 77 – Assay Info

List of PHASE COLOURS



Each Phase has a colour. Click on  icon (Fig. 76) to display which phase of the assay they represent (Fig. 78). This list only display phases that correspond to current session.




Fig. 78 – List for phase colors

5.2.2.2. Displaying Session Results

At the end of each session the results are automatically shown for immediate consultation. (Fig. 79).





Fig. 79- Results window

Select **Cancel** to close Results window without displaying results.

The **Error** button shows the following window in case the lid is opened during the reading:



Select **Detail** to open display with results for each analysed strip as shown here:




fig. 80– Details of Strip Results

Click icon to show print preview and further Detail for each strip (Fig. 81).




fig. 81– Display Details

In this display the fields concerning interpretation of results and additional notes (coloured Orange) can be changed before results are stored.

Click **Save** icon to save any changes – and Print preview displays how results can be printed and Print icon directly proceeds to print out results.

5.2.3. RUN LAST SESSION

The last run session can be run again from Main Menu, select **Session** and then **Run last session** from pull down list




Fig. 82– Run last session

This directly opens the Run Session display.

5.2.4. LOAD FROM HOST

This function concerns importing from Host and/or BCR (Fig. 83).




fig. 83

5.2.5. EXIT

From **Session** pull down menu select **Exit** to leave the programme.
If a session is running then a confirmation window appears first.




Fig. 84 – Exit

5.3. MENU DATABASE

5.3.1. ASSAY

See Paragraphs 6.2.

5.3.2. LAYOUTS

This display shows the different possible strip set-ups associated with different types of assays and is in two sections. (Fig. 55): The first (top) section Assay Type is for selecting Assays from a pull down listing. The selection shown here is Dot Blot. The second section (Layouts) displays a list of possible default set-ups associated with the Dot Blot assay and a series of options to access different functions for set-up (like for the assays).

Double click on any part of listing to display details of this set-up layout in a separate window and where the name of the Antigen can also be changed (Fig. 86).



APEX - [Layouts]

Assay type: Dot blot

Layouts:

Name	AntigenNumber
AD-GBM	1
AD-MID	1
AD-IFD	1
AD-ASCASD	1
AD-HND	2
AD-LIVER4	4
AD-CYT06D	6
AD-ENAD	6
AD-END	7
AD-FULLANAD	12
AD-LIVERD	5
AD-MPGD	3

Actions: New, Edit, Delete, Clone, Rename, Report, Protect, Import, Export

Detail: ☐ ☐ ☐

Buttons: Cancel, Save

fig. 85– Layout/Set-up




fig. 86 – Details of set-ups

5.3.2.1. New layout/Set-up

Click on **New** option in Layout display to create a new Strip set-up and select a number of antigens from the list displayed to coincide with the number of layouts (e.g.: AD-L01 has one antigen).





fig. 87– New set-up/layout

5.3.2.2. Changing Layouts

Click **Change** to change antigen name and controls, by clicking on **Chemical Data** tab, filling in **Antigen field**, display preview layout and possible print out by selecting **Report**.



The screenshot shows the 'APEX - [Layout detail]' window with the title 'AD-MID'. At the top, there is a 'Reading algorithm' dropdown menu set to 'ApdDbCrl'. Below this is a 'Well detail' section containing a horizontal bar with three circular indicators. To the right of this bar is a table with three columns: 'Position', 'Type', and 'Identifier'.

Position	Type	Identifier
1	Control	Positive
2	Antigen	M2
3	Control	Negative

Below the table is a horizontal scrollbar. Underneath the table is a tabbed interface with two tabs: 'Structural data' and 'Chemical data'. The 'Chemical data' tab is active, showing a 'Type' dropdown menu set to 'Control', a 'Score' input field with the value '100', and an 'Antigen' dropdown menu set to 'Positive'. At the bottom of the window are three buttons: 'Cancel', 'Report', and 'Save'.

fig. 88– Changing layout/set-up

5.3.2.3. Deleting a set-up/layout

Click on a layout from listing Fig. 85 and then select **Delete** option to delete this layout.

5.3.2.4. Cloning a layout

Click on a layout from listing Fig.90 and then select **Clone** option to copy this layout which is automatically renamed **NewLayout1** as seen here:



fig. 89 – Cloning layout

5.3.2.5. Renaming layouts

Click on a layout from listing (Fig.89) and then select **Rename** and then name can be changed directly on the listing.



fig. 90– Renaming layouts

5.3.2.1 Configuration protection

It is possible to protect a strip selecting it and pressing the button Protected.

A strip can be protected by the following user profiles:

- ☒ Admin
- ☒ Advanced Specialist

The remaining profiles cannot do it. The specialist can just create new configurations, import and export them, without the possibility to modify than ones created by the Advanced Specialist. The user can just to read the present configuration list.



Figura 1– Configuration protection

5.3.2.2 Report Print Out

Click on layout from listing (Fig.91) and **Report** to display a print preview of the selected layout.




fig. 91 – Report

5.3.2.6. Importing layouts

Click **Import** option to open a Windows display to load layouts in a *.xml file format.

5.3.2.7. Exporting layouts

Click **Export** option to open a display showing all current layouts. It is possible to de-select (un-click) Layouts not to be exported.

To export only a few layouts **Deselect All** And then re-select few layouts;
Select **OK** to open a Windows display to save layouts in a *.xml file format.




fig. 92 – Exporting layouts

5.3.3. RESULT ARCHIVES

To display previous analysis session results go to the Main Menu, select **Database** and then **RESULT ARCHIVES** from the pull-down menu (Fig.94)





Fig. 93 – Result archives menu

Both **Medical** and **Session reports** can be searched with filters: by Date and/or by Code for Medical Reports and only by Date to search through the Session Report

MEDICAL REPORTS



Patient ID	Date	Time	Session ID	Surname	Name
1	2009/05/20	09:39:28	New	a	a
2	2009/05/20	09:39:28	New	a	a

fig. 94 – Patient Search through medical reports

A date based search through the previous patient sessions is possible from the medical reports on this display. Input a date or period interval which may include the analysis sessions requested. Fill in the fields: **From** and/or **To** then click on **Search** to display a list of sessions processed during the period specified or from the date specified are then displayed.

To restrict the search further with an advanced search insert a sample code to find a previous patient analysis. A search can be started even by inserting only the Patient ID and/or Name and Surname then clicking **Search**.

The Search Results are displayed as shown below:

- ☐ PATIENT ID
- ☐ DATE
- ☐ TIME
- ☐ SESSION ID
- ☐ SAMPLE ID

SESSION REPORT

[illegible]

Fig. 95– Session report

Click on **Session Report** for a date-based search through previous patient sessions.

Input a date or period interval which may include the analysis sessions requested.

Fill in the fields: **From** and/or **To** then click on **Search** to display a list of sessions processed during the period specified or from the date specified are then displayed.

The Search Results are displayed as shown below:

- ☐ SESSION ID
- ☐ DATE
- ☐ TIME
- ☐ USER

5.4. MENU SETTINGS

5.4.1. LANGUAGE

The APEX software is configured for different language displays.

This option can be accessed from the **Settings** menu.

Select **Language** and then click on chosen language to be used in all displays as shown below (Fig.97)



Fig. 96 – Language

5.4.2. INSTRUMENT

Select **Settings** from the Main Menu and click on **Instrument** (Fig. 98) to access all information about the instruments connected.

This function is only active if licenses for more than one instrument have been loaded.




Fig. 97 – Instrument

5.4.3. USERS

A list of Users and their different software access levels can be created.




Fig. 98 – Users

From **Settings** in Main Menu select **Users** from pull down list to open **User List** as shown below (Fig. 100)




Fig. 99– User list

The **Admin** User can access all the software functions with no restrictions.
Select **New** to add any new Users.
Select **Modify** to change access and information of any User already listed.
Select **Delete** to cancel a User from the list.
Select **Change Password** to change the password of any User

5.4.3.1. NEW USER

To add a new User fill in the **New User Login Settings** as shown below (Fig.101)




Fig. 100 – New User Login

The following data is required:

- LOGIN Settings (obligatory)
 - ☐ User ID
 - ☐ User Password
 - ☐ Repeat User password
- USER CARD
 - ☐ User Registration number
 - ☐ User Name
 - ☐ User Surname
 - ☐ User role
 - ☐ Other Data
- ACCESS LEVELS (obligatory option boxes)
 - ☐ Simple User
 - ☐ Service
 - ☐ Specialist
 - ☐ Advanced Specialist

The programme manages four different pre-defined levels: **“SIMPLE USER”**, **“SERVICE”**, **“SPECIALIST”** and **“ADVANCED SPECIALIST”**.

Each User profile can access different functions:

When level has been selected (with a tick) click **SET**.

SPECIALIST:

SPECIALIST USER accesses all the same functions as the SIMPLE User and can also:

- **Set-up Assays (un protected)**

ADVANCED SPECIALIST:

The ADVANCED SPECIALIST can access all the same functions as a SIMPLE User and also:

- **Manage assays (protected and non)**
- **Manage Users**
- **Sets Cover Block level**



SERVICE:

SERVICE USER can access the same functions as the Simple User and:

- *Set-up instrument connections*

SIMPLE USER:

All the remaining functions can be accessed by Simple Users.

Please note: A single User can have more than one software Access Level to broaden possible functions accessed.

5.4.3.2. MODIFY

To change any data input for Users select User Name from User List.

The New User display appears where the LOGIN cannot be changed. After changing all the rest of the fields click on **SET** and the changes will be stored.

5.4.3.3. DELETE

To Delete a User, select it from list and click on Delete. A Message asking for confirmation appears before deletion.

5.4.3.4. CHANGE PASSWORD

To change a password select the USER then click on **Change password** to open display to input current password and new password. A Message asking for confirmation appears before password is changed.

5.4.4. CHANGE USER

To change User when software is already accessed and without relaunching it select Change User.




Fig. 101 – Settings - Change User

Select **Settings** from the Main Menu and click on **Change User** to access window as shown below. Type in New User login and password (fig. 103):



Fig. 102– Change User

Then Click **Enter**.

A confirmation window opens as shown in Fig. 104: click **Yes** if User changes level or **No**. Warning! Change from Admin User to User?




Fig. 103 – Confirm User change

Click **Cancel** (Fig. 103) or **No** from confirmation (fig. 103) to cancel all procedure.

5.4.5. CHANGE PASSWORD

Each User must have a personal password to access the software. In this way the password authenticates the User access levels and their different functions.




Fig. 104 – Settings - Change password

Here is the Change Password display. To change an existing password first input current password and then the new one which has to be input again for confirmation.

If the wrong confirmation password is inserted a Warning Message will appear to warn of error.



 A screenshot of the 'Change password' dialog box. It has a title bar with a close button. Inside, there is a 'User login' field with the text 'admin'. Below it is a 'Password' section containing three input fields: 'Current password', 'New password', and 'Repeat the password:'. At the bottom, there are two buttons: 'Cancel' and 'Set'.

fig. 105– Changing a password

Select **Set** to store password or select **Cancel** to avoid any effects of this function leaving old password active.

5.5. LABORATORY DATA

The APEX software can also include information about the laboratory using the instrument.

All information can be input by selecting **Laboratory Data** from Settings to open the **Laboratory data** window (Fig. 108). Any information will be printed out with all the analysis results.




Fig. 106– Laboratory data




Fig. 107 – Laboratory data

The details in this window are optional. Click Save option to store input data for possible display on reports.

5.5.1. DOCTOR DATA

The APEX programme also includes information about the Doctor conducting the tests: from **Settings** pulldown menu, click **Doctor Data** to open the **Doctor Data** (showing information about the Doctor) window (Fig. 110). Any information from this window is then printed out with all the analysis results.



Fig. 108 - Doctor data




Fig. 109 – Doctor data

Any information typed in this section is optional.

5.5.2. IMPORT SESSION DATA

To import any previously exported sessions exclusively in *.mdb file format.
Select "**Import Session Data**" and a display warns how importing new data could delete or overwrite previous sessions and asks confirmation to proceed.
Select "Yes" to open a Windows dialogue display to choose file.
A confirmation message appears when Import procedure is completed.




Fig. 110 – Import Session Data

5.5.3. EXPORT SESSIONS DATA

To export sessions only in *.mdb files, select **Export Sessions** and a Windows dialogue window opens to list possible destination paths where to save file.




fig. 111 – Export Sessions data

5.5.4. SETTING COMMUNICATION PORT

Open **Settings** from main menu and click on **Communication Port** to set up the serial port to use for communicating with the instrument (fig. 113).




fig. 112 – Settings – Communication port

The following **Communication port settings** opens (fig.114):




fig. 113– Setting Communication Port

Select Port from pull down listing as required (Com1, Com2 ...). In example fig.114 communication is for APBLOT instrument through Serial Port Com1.

Select **Set** to confirm settings as required; Select **Cancel** to stop procedure.

5.5.5. COVER BLOCK

This is only available for models with safety cover block functions and is optional.
Two different safety levels for blocking the cover can be set:

- ❑ High safety level: the session can only start when cover has been properly closed. It is blocked. Only if session is stopped or has reached the end can the cover be unblocked and opened.
- ❑ Low safety level: the session will start only after cover is closed but User can open the cover at any time during session (as in example shown in Fig. 116 below)

Cover Block movements (opening and closing) are recorded during sessions.
Security levels can be set-up for all different instruments.




Fig. 114 – Cover lock



Fig. 115 – Cover block settings

5.5.6. OPTIONS

From the pulldown menu **Settings** of the Main Menu click on **Options** (fig. 117).
(These options can only be accessed by Specialist and Advanced Specialist Users.)




Fig. 116 – Options

The **Options** menu can be used to set up different settings for all instruments and different parameters.


The **Options** display is divided into three: the **Instrument** section has a pull-down menu with a list of instruments.

This section is to set the same parameters for all instruments or one instrument at a time.

The **Parameters** section refers to the instrument(s) selected.

The **Detail** section shows the variables selected.

If all instruments are selected (as shown below) in the **Instrument** section, then the **Parameters** section displays the Type of error trace and set-up session form as in Fig. 118 below.



Options

Instrument: AllInstruments

Parameters:


Variable	Value
Type of error trace	LogFile
Setup session form	fmmWzdSesSetup

Detail:

Cancel Save

Fig. 117 – Options

The **Type of Error Trace** parameter is to display errors: only on screen (LogFile), with a message displayed (Msg) or both (MsgAndLogFile) or not displayed at all (None) (fig.119).




The screenshot shows the 'Options' dialog box with the 'Type of error trace' dropdown menu open. The menu options are: None, LogFile, Msg, and MsgAndLogFile. The 'MsgAndLogFile' option is currently selected.

Variable	Value
Type of error trace	MsgAndLogFile
Setup session form	frmWzdSesSetup

fig. 118– Type of Error Trace

Set Up Session Form (in detail section) can be used to display set up session window (assay-sample association) through Tree Graphics (frmSesSetup) or Wizard (frmWzdSesSetup).



The screenshot shows the 'Options' dialog box with the 'Setup session form' dropdown menu open. The menu options are: frmSesSetup, frmWzdSesSetup, and frmWzdSesSetup. The 'frmWzdSesSetup' option is currently selected.

Variable	Value
Type of error trace	MsgAndLogFile
Setup session form	frmWzdSesSetup

fig. 119– Set up session Form

Select Cancel option to exit Options display without making any changes. Select Instrument from pulldown listing and Parameters are displayed: Time arrangement Algorithm, Type selection Algorithm, Dilutor Final rinse time, Duration initial WASH prime, Duration first Buffer prime, Prime time after tank change, Duration of movement test prime Break on Error Mode if no liquid, Duration of Temporary interruption for Error (fig.121).




fig. 120– Options – Selecting single instrument

To choose different time **Time Arrangement Algorithm** is to set up different time intervals depending on requirements and session type, with assays that are 'similar' – which is when they have the same number of phases, same incubation time and same number of wash cycles.






fig. 121– Time arrangement algorithm

If assays are coherent then following selections can be made:

- CurrentSequence (this is for coherent and non-coherent assays): fitting the coherent assays together as if they were one, and leave non-coherent assays to the end.
- BlockShiftSingleDevice (non coherent assays): brings together assays where possible considering their different requirements.
- BlockPhaseShift (this is for coherent and non-coherent assays): this is to fit phases of coherent assays together (same phases).

Type Selection Algorithm can be used to change the time arrangement of the Assays (Kit number) as required: in automatic mode (**Automatic**), automatic and displaying a message (**AutomaticWithMessage**) or by interrupting a session (**BreakOnError**).

This parameter is only used if Current Sequences is also selected as "Time Arrangement Algorithm" whereby:

- **Automatic** is selected and finds non-coherent assays and treats them as **BlockShiftSingleDevice** without displaying any message
- **AutomaticWithMessage** can behave like Automatic parameter if requested.
- **BreakOnError** overrules **BlockShiftSingleDevice** when assays are non-coherent and the sessions are interrupted.




fig. 122– Type selection algorithm

Dilutor Final Rinse Time can set up how long a Prime lasts at the end of a session. This field can be edited.





fig. 123– Dilutor Final Rinse Time

Duration of initial Wash time is to set-up how long Prime lasts from Wash bottle at the beginning of a session. This field can be edited.




The screenshot shows the 'Options' dialog box with the 'Instrument' dropdown set to 'AP-Blot (s.n. 99)'. The 'Parameters' section contains a table with the following data:

Variable	Value
Time arrangement algorithm	CurrentSequence
Type selection algorithm	BreakOnError
Dilutor final rinse time (secs)	10
Duration of initial WASH prime (secs)	40
Duration first buffer prime (secs)	20
Prime time (secs) after tank change	20
Duration of movement test prime (...)	30
Break on error mode in case of mi...	WithTemporaryMessage

The 'Detail' section shows the 'Duration of initial WASH prime (secs)' field with the value 40. The 'Cancel' and 'Save' buttons are at the bottom.

fig. 124– Duration of initial WASH

Duration first buffer prime is to set up how long the prime lasts from buffer bottle at the beginning of a session. This field can be edited.



The screenshot shows the 'Options' dialog box with the 'Instrument' dropdown set to 'AP-Blot (s.n. 99)'. The 'Parameters' section contains a table with the following data:

Variable	Value
Time arrangement algorithm	CurrentSequence
Type selection algorithm	BreakOnError
Dilutor final rinse time (secs)	10
Duration of initial WASH prime (secs)	40
Duration first buffer prime (secs)	20
Prime time (secs) after tank change	20
Duration of movement test prime (...)	30
Break on error mode in case of mi...	WithTemporaryMessage

The 'Detail' section shows the 'Duration first buffer prime (secs)' field with the value 20. The 'Cancel' and 'Save' buttons are at the bottom.

fig. 125 – Duration first buffer prime

Prime time after tank change is the setting for choosing how long Buffer prime will last before each wash from this source.




fig. 126– Prime time after tank change

Duration of movement test prime is to choose how long prime lasts with WASH/Buffer depending on Prime Test set up from Tools menu.





fig. 127– Duration of movement test prime

Break On Error Mode is to set up instrument behaviour when there is no liquid - the absence of liquid can either be just indicated and the instrument continues procedure without interruptions (option is '**None**') or a message is displayed (**WithFixedMessage** option) with the choice of proceeding or interrupting the procedure or finally, when no liquid is found there is a pause time previously set up by the User (see later) after which the session continues with no further external interruptions (**WithTemporaryMessage** option).



Options

Instrument: AP-Blot (s.n. 99)

Parameters:

Variable	Value
Time arrangement algorithm	CurrentSequence
Type selection algorithm	BreakOnError
Dilutor final rinse time (secs)	10
Duration of initial WASH prime (se...	40
Duration first buffer prime (secs)	20
Prime time (secs) after tank change	20
Duration of movement test prime (...)	30
Break on error mode in case of mi...	WithTemporaryMessage

Detail:

Break on error mode in case of missing liquid: WithTemporaryMessage

Options: None, WithFixedMessage, WithTemporaryMessage

Buttons: Cancel, Save

fig. 128– Break on Error Mode

Duration Temporary Error Interruption is to set up how long an interval follows when an error in the absence of liquid is found when set up **WithTemporaryMessage**

5.6. TOOLS MENU

5.6.1. MOVEMENT TESTS

Select Tools from the Main Menu and then click on Movement Test (Fig. 130) from pull down list to access the following functions:

- A. Moving Probes in X,Y, Z directions and then moving the shaking plane
- B. Dilutor buffer
- C. Camera Optical test




Fig. 129 – Movement test

Open the **Movement Test** window (Fig. 131) and then select the instrument to be tested from pull down list. Select Start to start the test.

When test is completed the results are displayed next to each phase.




Fig. 130 – Movement test

5.6.2. PRIME TEST

Select **Tools** menu from Main Menu to access Prime Test function from pull down list (Fig. 132) and rinse hydraulic circuit.

From **Prime** Test window (Fig. 133) choose instrument to be primed and with which **Dilutor liquid** from pull down lists:

- Washer bottle (distilled water)
- Buffer n°1

Click **Start Prime** to begin. Duration of Prime is in seconds as set out by Duration of dilutor prime tool setting from the Options menu (see Fig. 133)




Fig. 131 – Prime test




Fig. 132 – Prime

5.6.3. LIGHT ON

The light can be turned on or off by the Specialist and Advanced Specialist User from the Main Menu Tools pull-down list (Fig.134).




Fig. 133 – Light On/Off

5.7. MENU ?

5.7.1. Guide

Select **?** from the Main Menu and the click on **Guide** from pull down list (Fig. 135) to display the Guide to the software in the language previously selected.




Fig. 134 – Guide

5.7.2. LICENSES

To manage Licenses select **?** from the Main Menu and click on **Licenses** from pull down list (Fig. 136).




Fig. 135 – Licenses

This opens Licenses window (Fig. 137) to add new licenses and display any current licenses if present.




Fig. 136 – Displaying license (s)

Right click on mouse on License as shown in Fig. 138 to delete license or display selected License expiry date.




Fig. 137 – Deleting/License displaying expiry date

5.7.3. INSTRUMENT INFORMATION

Open ? from the Main Menu and then click on **Instrument information** (Fig. 139) from pull down list to access all information about instruments connected (Fig. 140) including:

- Instrument Name
- Serial Number
- Instrument Code
- Devices
- Firmware Version




Fig. 138 – Instrument information




Fig. 139 – Instrument Information Report

5.7.4. SOFTWARE COMPONENTS

Select ? from the Main Menu and then click on **Software Components** from pull down list (Fig. 141) to access all information about installed software components (Fig. 142) which includes a list of files with each different version.




Fig. 140 – Software components




Fig. 141 – software components list

5.7.5. ABOUT

Select ? and then **ABOUT** from pull down list (Fig. 143) to display both manufacturer and software details (Fig. 144) including software version and manufacturer's contact details.




Fig. 142 – Information




Fig. 143 - Software and manufacturer contact details

SECTION 6 – SPECIALIST ACCESS

6.1. SPECIALIST and ADVANCED SPECIALIST

The Specialist and Advanced Specialist can access all Simple User functions.

The Specialist can also add unprotected assays to the assay list.

The Advanced Specialist can also add both password protected and unprotected assays, change users and also access and change cover block security levels.

6.2. SETTING UP A NEW ASSAY

To set up a new assay select **Database** from the main menu and click **Assays** from pull down list (Fig.145 below).




Fig. 144 – Assays





Fig. 145 – Assays

This opens the **Assays** window divided into two parts which show:

- All the stored assays with relevant information
- A series of options (on the right side) for different functions (described later).

The following displays are examples of what kind of information needs to be input to set up a Dot Blot assay.



The screenshot shows a 'New Assay' window with a dark background and a grid of wells. The window has tabs for 'DATA', 'KIT', and 'RESULT'. The 'DATA' tab is active. The following fields are visible:

- Assay name: A2
- Short name: A2
- Assay group: G1 (dropdown menu)
- Layout: AD-ASCASD (dropdown menu)
- Washing buffer: Buffer
- Protection: ☐
- Assay type: DotBlot
- Created: 15/07/2009 17:00:32
- Last change: 15/07/2009 17:00:32

A 'Save' button is located at the bottom right of the window.

Fig. 146 – Information for Dot Blot Assay

- Click **New** option  to set up a new assay from **New Assay** window and then fill in relevant fields:
 - ❑ Select **Assay group** from pull down listing which shows all the possible groups.
 - ❑ Type in **Assay Name** and **Short Name** to assign to assay for easier identification (A1 in the example).
 - ❑ The **Wash Buffer** field has a default value that is the name of the buffer required for washings; such name can be changed as required.
 - ❑ The **Protection** box option can be selected to set up limited access to this Assay Layout for only Advanced Specialist level Users.
 - ❑ The **Assay Type**, **Last Change** are already compiled
 - ❑ Select type of Strip **Layout** from pull-down listing to associate with assay or click Browse  . Option to open APEX Layout Display as shown in Fig. 18 and choose from list of possible layouts.




fig. 147– Layouts

- Click on Report  to display selected strip layout for the assays (Fig 148).





Fig. 148 – Layout

- Select SET option  to store Layout for this New Assay.

6.2.1. Phases Settings

After the **DATA** display has been filled in for the New Assay set-up, click on KIT Tab along top left of New Assay window to open the KIT information for the Assay parameter settings display.



Controls

Dispensation mode: Begin

Number	Value	Vial

Value:

Vial: LD14H30

Work Session

Phase:

Vial:

Number	Phase	Vial
1	Pretreat	D34H67
2	Sample	
3	Conjugate	D34H67
4	Substrate	D34H67
5	Reading	

Save

fig. 149– Setting up Assay Kit

The KIT display is in two sections: the top section shows information about **Controls** and the bottom section shows the Phases and different Reagents for **Work Sessions** or Phases.

ADDING A NEW CONTROL

The **Controls** section shows the following information in a table:

- ☐ Control Number (increasing numbers)
- ☐ Value
- ☐ Vial

Only the Value can be selected from the pull down list as the Control Vial is already pre-set in the only position in the Reagent Rack.

The Values can be:

- Positive
- Negative
- Cut Off

Select **New Control** from pull down list to edit the **Value** section shown below and input new information.

The screenshot shows the 'New Assay' window with three tabs: DATA, KIT, and RESULT. The 'Controls' section is active, showing a 'Rack' dropdown set to 'None' and a 'Dispensation mode' dropdown set to 'Begin'. Below these is a table with 3 rows and 3 columns: Number, Value, and Vial.

Number	Value	Vial
1	Positive	CD14H38
2	Negative	CD14H38
3	Cut-Off	CD14H38

Below the table are 'Value' and 'Vial' dropdown menus, both set to 'CD14H38'. There are 'Add' (left arrow) and 'Delete' (X) buttons. The 'Working Session' section at the bottom has a 'Phase' dropdown set to 'D34H67' and a 'Vial' dropdown set to 'D34H67'. It also has 'Add' and 'Delete' buttons. A table with 5 rows and 3 columns (Number, Phase, Vial) is shown below:

Number	Phase	Vial
1	Pretreat	D34H67
2	Sample	
3	Conjugate	D34H67
4	Substrate	D34H67
5	Reading	

Navigation arrows (up, down) are next to the table. A 'Save' button is at the bottom right.

Fig. 150 – Setting up Kits

Select **Add** option to store set-up which is immediately displayed on Controls table;

Select **Delete** option to cancel selected control from Controls table.

ADDING/DELETING A PHASE

In the Work Session section the table shows:

- ☐ Number (increasing number)
- ☐ Phase
- ☐ Vial

A Phase can be selected from Phase pull down listing and currently only one type of Vial can be selected.

Select **Add** option to store information and display them in the table;

Select **Delete** option to delete previously selected phase from table.


The options move up and down the order of Phases and can be used to change the chronological order of the phases to be run.

Select a Phase from the table and click **Down** arrow and the Phase will move down to the following position.

Select a Phase from the table and click the **Up** arrow to move the Phase up to the previous position

CHANGING PHASE SETTINGS

Double click one of the Phases in the Phase table (Fig. 152) to open the following Sub-Phase window for each Phase.



Pretreatment

Thermoregulation Temperature (°C)

☒ **Dispensation**

- Dispensed volume (µl):
- Initial aspiration:
- Dispensation mode:
- Probe washing volume (µl):
- Aspiration speed (Hz):
- Dispensation speed (Hz):
- Waste volume (%):
- Bubble volume (µl):
- Shaker state:
- Shaker speed (Hz):
- Shaker slope:
- Aspiration slope:

☒ **Wash**

- Dispensed volume (µl):
- Washing cycle number:
- Soak time (s):
- Source:
- Wash mode:
- Dispensation speed (Hz):
- Shaker state:
- Shaker speed (Hz):
- Shaker slope:
- Aspiration slope:

☒ **Incubation**

- Time interval (min):
- Shaker state:
- Shaker speed (Hz):
- Shaker slope:

☒ **Cleaning solution cycle**

- Aspiration speed (Hz):
- Dispensation speed (Hz):
- Volume (µl):

fig. 151– Sub-Phase

Each Sub Phase is made up of three stages in order:

- ☐ DISPENSATION
- ☐ INCUBATION
- ☐ WASH

The Stop is the only default phase which has un-selected the incubation and wash phases because they cannot be set-up.

WASHING PROBES BY SOLUTION shows parameters for washing the probes.

All other sub-phase parameters can be set-up by filling in the text fields or selecting from pull-down lists.

The following parameters are required in the DISPENSATION section of the display:

- Dispensed volume (µL)
- Initial aspiration (µL)
- Dispensation mode
- Probe washing volume (Hz)
- Aspiration speed
- Dispensation Speed (Hz)
- Waste volume (µl)
- Bubble Volume
- Shaker State
- Shaker Speed (Hz)
- Shaking Slope
- Aspiration slope

The following parameters are required for the INCUBATION section:

- Time interval (min)
- Shaker State
- Shaker Speed (Hz)
- Shaking slope

The following parameters are required in the WASH section:

- Dispensed volume (µl)
- Washing cycle number
- Soak time (s)
- Source
- Wash Mode
- Dispensation speed (Hz)
- Shaker State
- Shaker speed (Hz)
- Shaker Slope
- Aspiration slope

☐ The following parameters are required in the Cleaning solution cycle

- Aspiration Speed (Hz)
- Dispensation Speed (Hz)
- Volume (µl)

A box is present next to each SubPhase: if the box is ticked ☒ then the relevant information has been inserted. Un-tick the box to delete the SubPhase.

To access the **Controls** and **Samples** as shown in Fig. 153 below, the **Predilution** section has to be filled in with the following information:

- Dilution Ratio
- Diluent
- Type
- Probe Washing Volume (µl)
- Diluent aspiration speed (Hz)
- Diluent dispensation speed (Hz)
- Liquid aspiration speed (Hz)
- Diluent bubble volume (µl)

To access the parameters of the other SubPhases then click on the relevant tab (Incubation/Wash) along the top of the display.





Fig. 152 –Sample and Controls phases

If at least one control has been inserted in the Samples and Controls display then a tick will appear ☒ for the Controls so that dispensation and predilution settings can be filled in.
Please note: Sample predilution and Controls settings are independent.



The screenshot shows the 'Samples and Controls' window with the 'Predilution' subphase selected. The interface is divided into two main sections: 'Dispensation' and 'Predilution'. Each section has a 'Samples' and 'Controls' column with checkboxes and input fields.

Parameter	Samples	Controls
Dispensed volume (µl)	1500	1500
Previous aspiration	Yes	Yes
Dispensation mode	1	1
Probe washing volume (µl)	2000	2000
Aspiration speed (Hz)	480	480
Dispensation speed (Hz)	1633	1633
Waste volume (%)	10	10
Bubble volume (µl)	12	12
Shaker state	On	On
Shaker speed (Hz)	1400	1400
Shaker slope	6	6
Aspiration slope	6	6

Parameter	Samples	Controls
Dilution ratio	1/100	
Diluent	Diluent1	
Type	On fly	
Probe washing volume (µl)	2000	
Diluent aspiration speed (Hz)	930	
Diluent dispensation speed (Hz)	1633	
Liquid aspiration speed (Hz)	480	
Diluent bubble volume (µl)	12	

fig. 153– Predilution SubPhase – samples and controls

For the **Reading** Phase only input amount of seconds for Minimum Delay before Reading setting.




fig. 154– Reading Phase

N.B.: If the specific assay set up does not include any of these Sub-Phases then the relevant parameter boxes should be disabled by deleting check.

All values and information input are saved when each Phase display is closed (click Exit)

RESULT

Select **Result** from the tabs along the top of the Assays display to access the options to choose the different classes for presenting the Results (fig.156).

Each class defines a range of values to diagnose and evaluate data. A lower and upper limit are defined for each class, where the upper limit of the one class is defined by the lower limit of the previous class.




Fig. 155 – Results

ADDING/DELETING A NEW CLASS

In the Result area the following table appears to input:

- ☐ Interval
- ☐ Lower limit
- ☐ Upper limit

To set up a new class click on NEW  compile the **Interval** text box and **Lower limit** fields and the click on **Add**.



All selections are set and displayed on the table depending on the Lower Limit Field.

Select **Cancel**  to delete previous selected phase in the table.

Select **Save**  to store settings – the new assay will be now displayed in the **Assays** list (Fig. 157).




fig. 156– New Assay displayed

6.2.2. Changing an Assay

From the **Assay** display (fig. 158) select the assay to be changed and then select **Edit** (fig. 159).



fig. 157– Assays

Or click two times on name from list.




Fig. 158 – Changing an assay

Input different parameters in the Modify Assay display and then click on **Save** to store new information.

Deleting an assay

To delete an assay: select the assay from the database list as shown in Fig. 158 and then click **Delete**



A confirmation window appears before actually Deleting an assay: click Yes to definitively cancel or delete the assay from the list or click NO and the assay remains stored.

6.3. Copying an existing assay

To copy an assay already stored – select this assay from the database list as shown in Fig. 158 and then click



Copy.

The newly 'generated' Assay will appear in the assay list with the name of the copied assay followed by an automatically generated number.

To change this name select this Assay from the list, then click Edit to insert the new name and new short name.

6.4. Printing out an assay

To print out assay parameters, select the assay from the database list which appears in Fig. 158 and then click



Print

A normal print preview Windows dialogue window appears showing the normal range of print options.


6.5. Importing a protocol

Any previously stored Assay exported in a *.xml format can be Imported.

All assays in the Protocol to be imported will be added to Assays already stored.

If name of imported protocol is the same as an existent assay protocol then a message will appear asking whether old assay should be substituted with the new assay or not.

If no substitution is desired then re-name original Assay and then proceed with Importing Assay.



From Assay display (Fig. 128) select **Import** to open possible selection to import.

A typical Open File window appears and the selected assay destination file path can be input as a *.xml type file.




fig. 159 – Import Dialogue Windows display

When desired protocol is defined then select **Open** to open the **Assay Import/Export** window where one or more assays can be selected for import or export.

Click next to Assay name to select Assays to be imported or exported and then **OK**.

Click **Deselect all** to cancel Import/export procedure.




fig. 160- Dialogue display

6.6. Exporting an assay

Assays can also be saved in an *.xml file to then store and/or import them later.

Click **Export**  (see Fig. 162) to open the **Assay Import/Export** window where one or more assays can be selected for import or export.




fig. 161– Import/Export assay window

After selecting Assays click **OK** and a window will appear to type in the destination address of where the file or files should be exported to and stored.

The exported assay is saved as an *.xml file which automatically includes all the information about the stored assay.

Click  to de-select with one click all assays in this protocol, and then choose the few to be exported.

SECTION 7 – MAINTENANCE



Only Qualified personnel are entitled to carry out maintenance (see Section 1 of this manual). Any maintenance should follow general rules and safety precautions described in Section 1 of this manual.



The instrument is used with chemical reagents and other dangerous CHEMICAL SUBSTANCES (corrosive, irritant and harmful) which can damage health. When this label is present particular care and attention must be taken to follow all manufacturer recommendations.



This instrument handles potentially infected samples (urine, human serum etc.). INFECTIONS or CONTAMINATION might occur. Always follow general safety precautions in the presence of such biological substances. Use protective clothing, gloves and goggles.

7.1. MAINTENANCE CHECKS & PRECAUTIONS

The suggested maintenance measures are indicated in this chart:

Type of service	When
Cleaning	Daily and / or after a work session
Inspection and checking all parts	Recommended period between checks is three to six months depending on use.
Movement Tests	Every 6 months

To maintain instrument in optimum working conditions regular servicing should be carried out **as indicated in the table** above. The instrument should also be serviced after a long period of non-use. After particularly heavy use the maintenance checks and procedures described above should be carried out more regularly.

7.1.1. Cleaning



Warning: when cleaning do not use alcohol or similar solutions!

While cleaning do not wet any electrical parts (connectors etc.). If some electrical parts accidentally get wet, carefully dry these parts before turning on and using instrument.

Thorough cleaning of the instrument is recommended after each work session:

- Follow decontamination procedures described in Appendix A or other equivalent internal standard procedure
- To clean exterior of instrument use a water-damp cloth with a little detergent
- Keep all instrument surface work areas clear and clean by removing any dirt or possible drips.

7.1.2. Instrument testing

Regular instrument testing is required as follows:

- Check functions of control panel
- Line testing

7.1.3. Inspections and checks

All instrument parts should be inspected and checked regularly - for instance:

- Inspect power cable, substitute if damaged/deteriorated
- Check hydraulic dispensation circuit and substitute any connecting tubing if necessary

7.2. REPAIRS

User(s) should ensure that only regular ordinary maintenance is carried out. If any faults or problems occur while working with the machine, only request a Qualified Technician to service the instrument. **REPAIRS (like additional maintenance) see Section 1 of this manual – should be carried out by:**

- ☐ **Only Qualified Technicians** from Authorised Centers (specifically trained in this specific instrument servicing) should carry out any REPAIRS and instrument service maintenance, with specific Service Technical Manual and original spare parts.



Alterations or modifications of the instrument are not allowed. The owner/user is liable for any improper use or modification of the instrument, in addition to any consequences.

SECTION 8 - PACKING, TRANSPORTATION, STORAGE AND PUTTING INSTRUMENT “OUT OF SERVICE”

8.1. GENERAL INFORMATION

For packing, transport and storage procedures the instrument has to be put “out of service”. Carefully follow instructions below:



The instrument has been used with potentially infected samples (urine, human serum etc.). **INFECTIONS** and **CONTAMINATION** can occur.

General safety precautions should be carefully observed when dealing with potentially infectious biological substances.

Before putting the instrument out of service **IT MUST BE DECONTAMNIATED!**

(See APPENDIX A).

Before instrument transport or storage, complete **DECONTAMINATION DECLARATION** with correct date and signature by a qualified USER.

8.2. PUTTING INSTRUMENT “OUT OF SERVICE”/DECOMMISSIONING

Remove instrument power cable from PC and device.

If the instrument has been used - **disinfect and clean properly** - following disinfection instructions and Internal Laboratory procedures.



If the instrument has to be moved or stored, enclose and/or display **Decontamination Certificate** (see Appendix A) dated and signed by a qualified User.


8.2.1. Temporary Storage

Deactivate instrument and use specific warning stickers for “OUT OF SERVICE MATERIAL” on storage container

8.3. PACKING THE INSTRUMENT

Careful instrument packing is required when any transportation/delivery is involved. To pack the instrument follow instructions described below:

- Position instrument inside the original case.
- Apply proper stickers and safety indications on the box:



8.4. TRANSPORTATION

Manually transporting the instrument without its case should be limited to the work area, avoiding unnecessary knocks or collisions.
Original packing can be used for other transport/delivery.



During transportation the instrument must be **accompanied by the Decontamination Certificate** (see Appendix A) preferably visible on the outside of the original packing.

8.5. STORING THE INSTRUMENT

For instrument storage for any period of dis-use, use original packaging if possible as mentioned earlier.
The storage area should not be damp (80% max. humidity level) and temperature range between 4-60 °C.



The instrument has to be accompanied by **Decontamination Certificate** (see Appendix A) possibly displayed clearly on storage container.



Safety warnings and general rules have to be observed when moving or lifting the instrument.

8.6. DISPOSAL PROCEDURES



The manufacturer complies with EU directives: 2002/95/CE, 2002/96/CE and 2003/108/CE. These regulate how to reduce using dangerous substances present in electrical and electronic devices and their disposal. The above symbol (cross over dustbin) either on the instrument or on the packaging shows that after end use the product must be disposed of separately and in compliance with regulations about these products.

The manufacturer of this instrument is responsible for disposing of it according to the principles of differentiated waste disposal.
Users should contact the manufacturer and follow indications to enable its disposal at the end of its usable life.

This kind of approach: whether for re-cycling purposes, treatment and environmentally sustainable disposal contributes to avoiding possible negative effects on the environment and encourages re-use and/or recycling of material used for instrument manufacture all or in part.
Any illegal disposal of this instrument can result in administrative penalties as set out by current legislation and regulations.



The instrument has to be accompanied by **Decontamination Certificate** (see Appendix A) displayed clearly on storage container.

SECTION 9 – APPENDIX

APPENDIX A

A.1. INSTRUMENT DECONTAMINATION



The instrument involves the handling of potentially infected samples (urine, human serum etc.)
INFECTIONS or CONTAMINATION might occur.
Every part and accessory of the instrument must be considered potentially infected.

A.1.1. DECONTAMINATION PROCEDURE

Use one of the following wide band disinfectant solutions (or equivalent):

- **Lysetol**
- **Aseptisol**
- **Germocid Plus**



Only use solution for instrument surfaces.

1. Wear disposable gloves, protective goggles and suitable clothing
2. Prepare a sealable bag for disposable items used for the decontamination and label and seal bag accordingly
3. Unplug instrument from electrical sources to avoid explosions
4. Remove all accessories and disinfect those which are included and packed with the instrument
5. Spray disinfectant solution on instrument surfaces or use a cloth or paper soaked in a disinfectant solution to wipe instrument surfaces
6. Leave solution on instrument for 10 minutes and repeat step (5)
7. Leave solution on instrument for 5 hours, clean off instrument surface with a light detergent or distilled water to eliminate any dirt or disinfectant solution
8. Carefully dry instrument
9. Put instrument and accessories into original box/packaging
10. Wash and disinfect the hands with a light detergent
11. Fill in the Decontamination Declaration and enclose with instrument

A.2.1. Decontamination Declaration

The Decontamination declaration (shown below) IS OBLIGATORY – correctly compiled and enclosed with the instrument before storage or shipping for maintenance or repair service.
The declaration has to be VISIBLY displayed on instrument packaging.



DECONTAMINATION DECLARATION

Instrument: _____

Model: _____

S/N: _____

We declare that:

- ☐ instrument and accessories have never come into contact with dangerous biological substances,
 - ☐ or
- ☐ instrument and accessories have been decontaminated to disable or eliminate any potentially dangerous biological substance for personnel.

Client name: _____

User: _____

Address: _____

Date _____ Signature

APPENDIX B

B.1. HYDRAULIC CIRCUIT DIAGRAM

